



## BriteMED® pCART-N1


Ergonomic hospital trolley for medication purpose with front drawers

### Features

- Aluminium-alloy structure ensures good durability and bump-resistant during movement
- 3-layer medication drawers allows multiple medicine storage and classification
- Powder coated surface facilitates clean and sterilization purpose
- Optional VESA 75/100 arm enables quick adjustment for ideal viewing position
- Concave desktop design and scratch-resistance surface increase the durability

### Mobile Nursing Workstation

By integrating BriteMED® Battery Medical Cart PC, the computerized pCART-N1 becomes an intelligent medication workstation which creates clinical mobility that streamline the hospital workflow and enhance ward round efficiency.



Battery Medical Cart PC  
(POCm-W22C-ULT3)  
(POCm-W24C-ULT3)

VESA 75/100 Arm  
(7C000-LA6A2)

### Specifications


Model		pCART-N1
Description		Ergonomic for medication use with front drawers.
Size (Lx W x H)		68 x 54 x 103 cm (±1cm)
Structure		Aluminum-alloy
Weight		66.6 kg
Working Table ( ABS surface)	Dimensions (L x W)	59 x 45 cm
Side Table ( ABS surface)	Loading weight (max.)	8 kg
*Supported Arm	Loading weight (max.)	12 kg
Standard Drawer Configuration: * Drawers are metallic and powder coated, including stainless steel slippery track and drawer lock.	Total Height	27"
	3" keyboard drawer	1
	3" drawer	1
	3" medication drawer	3
	6" drawer	2
Wheels (4 dust-proof wheels with 2 in brake function)	Size	5"

\* Supported Arm is optional. Please see "Options" section for more information.

## Ordering Information

Part No.	Description
pCART-N1-N-W-A-U-R2	Mobile Cart N1 white color with ABS TOP, 3" keyboard drawer x 1, 3" drawer x 1, 3" medicine box drawer x 3, 6" drawer x 2, pill box x 15
pCART-N1-A-W-A-U-R2	Mobile Cart N1 white color with ABS TOP, 3" keyboard drawer x 1, 3" drawer x 1, 3" medicine box drawer x 3, 6" drawer x 2, pill box x 15, VESA 75/100 Monitor / Panel PC Arm x 1

## Options

	Part No.	Description
	PCA-VBH-R1	Bottle holder(31x11x8cm) with side rail
	PCA-VMS-R1	Multi-storage holder(33x13x18cm) with side rail
	PCA-VGD-R1	Glove dispenser holder(26x8x13cm) with side rail
	PCA-VCH-R1	Catheter holder (10x10x48cm) with side rail
	PCA-VIV-R1	IV pole
	PCA-VSC-3L-R1	Sharps container(yellow 3L x1) with holder(18x12x11cm) & side rail
	7C000-LA6A2	VESA 75/100 Monitor / Panel PC Arm

	Part No.	Description
	PCA-VSC-R1	Sharps container (red 1.7L x1) with holder (12x12x9cm) & side rail
	PCA-VSB-R1	Metal waste bin(10L x2) (40x17x34cm)
	PCA-VWC-R1	Plastic waste bin with lid(9L x1) & holder(26x15x18cm) & side rail
	PCA-DPFD3-R1	3" fixed divider set (short x3 & fixed x3 & strip x2)
	PCA-DPFD6-R1	6" fixed divider set (short x2 & fixed x2 & strip x2)
	PCA-DMB-R1	3" medicine bin (232x113x75mm)
	POCm-W22C-ULT3/ POCm-W24C-ULT3	22" & 24" Battery Medical Cart PC



pCART-N1-W\_  
datasheet\_V1\_20210902

### Britemed Technology Inc.

3F., No.306/306-3, Sec. 1, Datong Rd., Xizhi Dist., New Taipei City 22146, Taiwan  
 TEL: +886-2-86919498  
 FAX: +886-2-86919468  
[www.britemed.com.tw](http://www.britemed.com.tw)

Specifications may change without further notice.  
 © 2021 BriteMED Technology Inc. All rights reserved.