

BriteMED® pCART-N3

Compact design for documentation purpose with front drawers



Features

- Locked drawer ensures the safety of all properties
- Powder coated surface facilitates clean and sterilization purpose
- Optional VESA 75/100 arm enables quick adjustment for ideal viewing position
- Enhanced strength with aluminum-alloy structure and electro-welding metal sheet

Specifications

Model		pCART-N3
Description		Compact for documentation use with front drawers.
Size (L x W x H)		51 x 46 x 94cm (±1cm)
Structure		Aluminum-alloy
Working Table (ABS surface)	Dimensions (L x W)	45 x 39cm
Side Table (ABS surface)	Loading weight (max.)	8 kg
Supported Arm (VESA 75/100)	Loading weight (max.)	12kg
Standard Drawer Configuration: * Drawers are metallic and powder coated, including stainless steel slippery track and drawer lock.	3" keyboard drawer	1
	3" drawer	2
	3" medication drawer	N/A
	6" drawer	1
	9" drawer	1
Customized Drawer Selection		Total height 27", freely combining 3",6",9" drawers
Wheels (4 dust-proof wheels with 2 in brake function)	Size	4"

Ordering Information

Part No.	Description
PCART-N3-PP1-MC-R1	Medication Cart N3 pink color with 3" keyboard drawer x 1, 3" drawer x 2, 6" drawer x 1, 9" drawer x 1

Options

	Part No.	Description		Part No.	Description
	PCA-VBH-R1	Bottle holder(31x11x8cm) with side rail		PCA-VSC-R1	Sharps container (red 1.7L x1) with holder (12x12x9cm) & side rail
	PCA-VMS-R1	Multi-storage holder(33x13x18cm) with side rail		PCA-VSB-R1	Metal waste bin(10L x2) (40x17x34cm)
	PCA-VGD-R1	Glove dispenser holder(26x8x13cm) with side rail		PCA-VWC-R1	Plastic waste bin with lid(9L x1) & holder(26x15x18cm) & side rail
	PCA-VCH-R1	Catheter holder (10x10x48cm) with side rail		PCA-DPFD3-R1	3" fixed divider set (short x3 & fixed x3 & strip x2)
	PCA-VIV-R1	IV pole		PCA-DPFD6-R1	6" fixed divider set (short x2 & fixed x2 & strip x2)
	PCA-VSC-3L-R1	Sharps container(yellow 3L x1) with holder(18x12x11cm) & side rail		PCA-DMB-R1	3" medicine bin (232x113x75mm)
	7C000-LA6A2	Monitor/ Panel PC Arm			



BriteMED Technology Inc.

3F., No.306/306-3, Sec. 1, Datong Rd., Xizhi Dist., New Taipei City 22146, Taiwan

TEL: +886-2-86919498

FAX: +886-2-86919468

www.briteMED.com.tw

Specifications may change without further notice.

© 2020 BriteMED Technology Inc. All rights reserved.

Celebrate Humanity With Technology