

BriteMED® pCART-D1

Slim and compact medical trolley for documentation purpose

Features

- Dust-proof wheels and easy-to-clean surface
- Enhanced usability by included basket and device holder
- Easy adjustment of table height by air-pressure spring design
- Flexible to integrate panel PC and battery for e-medical solution
- VESA 75/100 arm enables quick adjustment for ideal viewing position




Specifications

Model		pCART-D1
Description		Slim and economical for documentation use in limited space.
Size (L x W x H)		52 x 54 x 150cm (±1cm)
Structure		Aluminum-alloy
Working Table (ABS surface)	Dimensions (L x W)	54 x 45cm
Supported Arm (VESA 75/100)	Loading weight (max.)	10kg
Standard Drawer Configuration: * Drawers are metallic and powder coated, including stainless steel slippery track and drawer lock.	3" keyboard drawer	1
Mouse Holder		1
Document Basket		1
Barcode Reader Holder		1
Wheels (4 dust-proof wheels with 2 in brake function)	Size	4"

Ordering Information

Part No.	Description
PCART-D1-W00-N-R1	Medication Cart D1 white color with height adjustable desktop 90~120cm, keyboard tray, barcode reader holder, document basket, mouse holder, 4" casters, VESA 75/100 arm


Flexible and Solid Computer Arm



Height Adjustable Table




Smooth Movement & Wheels with Brake Plate



Easy to Remove



Options

	Part No.	Description		Part No.	Description
	PCA-HCA-001-R1	3" height drawer with lock, 27x44cm		PCA-HCA-002-R1	Desktop PC holder with belt, 25x15x60cm

