

## BriteMED® pCART-D2

Slim-type documentation cart



### Features

- Easy-to-clean surface
- Slim design for space-saving advantage
- Dust-proof and user-friendly caster design
- Easy height adjustment for ideal viewing position
- VESA 75/100 arm for different mounting options


### Specifications


Model		pCART-D2
Description		Slim and economical for space-saving advantage
Size (L x W x H)		61.1 x 47.5 x 126.5cm
Working Table ( ABS surface)	Dimensions (L x W)	51.7 x 26.1cm
Supported Arm (VESA 75/100)	Loading weight (max.)	12kg
Standard Drawer Configuration: * Drawers are constructed from metallic foothold and plastic material, featuring stainless steel sliding tracks and a drawer lock.	3" keyboard drawer	1
Mouse Holder		1
Document Basket		1
Barcode Reader Holder		1
Wheels ( 4 dust-proof wheels with 2 in brake function)	Size	4"

### Ordering Information

Part No.	Description
PCART-D2-W00-R1	Medication CartD2 white color with height adjustable desktop 80~130cm, keyboard tray, mouse holder, 4 inch medical grade casters, VESA 75/100 arm.

### Options

	Part No.	Description
	PCA-HCA-004-R1	Medicine Container With Lock 42.6*41.4*10.8 cm

	Part No.	Description
	PCA-HCA-003-R1	Drawer With Lock 42.6*41.5*10.8 cm