

BriteMED® PCA-AUW10 Articulated Arm for Display Mounting




Features

- Support up to 12kgs
- Flexible to mount on wall
- Height adjustment from -30°(down) ~ +65°(up)
- Multi-angle adjustment by 180 degrees

Specifications

Structure		Steel/Aluminum Alloy/Plastic
Cable Management		Yes
Net Weight		1.8 kg (3.97lbs)
Dimensions (W*H*D)		410 x 370.5 x 111mm (16.14" x 14.59" x 4.37")
Hinge Module	Swivel	+/-90°
	Tilt	+/-20°
	Pivot	360°
Display Support	Display Type	Single Monitor / Panel PC
	Support Max. Weight	12kg (26.46lb)
Height Adjustment (Arm)	Height Adjustment Range	315 mm
	Tilt	-30°(down) ~ +65°(up)

Mechanical Drawing



Britemed Technology Inc.

3F., No.306/306-3, Sec. 1, Datong Rd., Xizhi Dist., New Taipei City 22146, Taiwan

TEL: +886-2-86919498

FAX: +886-2-86919468

www.britemed.com.tw

Specifications may change without further notice.

© 2024 BriteMED Technology Inc. All rights reserved.